

PentaFerte

MEDICAL DEVICES SOLUTIONS



PENTAFLUSH®

About PentaFlush® pre-filled saline syringes

PentaFlush® offer a complete range of sterile pre-filled syringes for the flushing of vascular access devices.

General specifications

- ✓ Sterile, single use medical device
- ✓ Suitable for sterile field
- ✓ Standard 10 ml barrel for 3, 5, 10 sizes to reduce internal pressure and avoid catheter damage
- ✓ **Class II B** medical device
- ✓ 3 years shelf life
- ✓ During the syringe manufacturing process the following substances are not used or intentionally incorporated into their materials: DEHP, PVC, LATEX, PRESERVATIVES

Components specifications

- ✓ **Barrel:** transparent polypropylene with plunger retention ring
- ✓ **Plunger:** opaque polypropylene
- ✓ **Plunger stopper:** synthetic rubber with two sealing rings
- ✓ **Lubricant:** medical grade silicone oil
- ✓ **Graduation scale:** 1ml interval
- ✓ **Saline solution:** sterile saline solution NaCl 0,9%
- ✓ **Unit packaging (sterile barrier system):** Medical Grade Paper & Plastic Film

Sterilization method

Gamma radiations, cycle validated according to EN ISO 11137.

Code	Description	Units/Carton	Boxes/Carton
060022424	3 ml pre-filled syringe 0,9% NaCl sterile	180	6 x 30 PCS
060022524	5 ml pre-filled syringe 0,9% NaCl sterile	180	6 x 30 PCS
060022624	10 ml pre-filled syringe 0,9% NaCl sterile	180	6 x 30 PCS
060022724	20 ml pre-filled syringe 0,9% NaCl sterile	80	4 x 20 PCS